



### **Product description**

Order number	65.1121
Product description	Screw cap, natural, suitable for tubes Ø 11.5 mm, 1,000 piece(s)/bag

### **Product characteristics**

Cap	standard screw cap
Closure type	screw cap
Suitable for	tubes Ø 11.5 mm

### **Size**

Diameter	13.5 mm
Height of product	8.7 mm

### **Material & colours**

Product material	High Density Polyethylene (HD-PE)
Colour of product	natural

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	390000
Depth of case	418 mm
Width of case	278 mm
Height of case	130 mm
Case volume	0.0151 cbm
Weight of product	0.0003 kg
Weight of case	1.91 kg
EAN of inner box	4038917085397
EAN case	4038917004077

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

