



### **Product description**

Order number	65.729
Product description	Screw cap, natural, suitable for tubes Ø 15.3 mm, 1,000 piece(s)/bag

### **Product characteristics**

Cap	standard screw cap
Closure type	screw cap
Suitable for	tubes Ø 15.3 mm

### **Material & colours**

Product material	High Density Polyethylene (HD-PE)
Colour of product	natural

### **Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Latex-free	yes
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	280000
Depth of case	428 mm
Width of case	383 mm
Height of case	175 mm
Case volume	0.0287 cbm
Weight of product	0.001 kg
Weight of case	5.45 kg
EAN of inner box	4038917083775
EAN case	4038917365239

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

