



**Product Specification**  
**Push cap, natural, suitable for tubes Ø 10 and 11 mm**

**Page 1**



**Product description**

Order number	65.782
Product description	Push cap, natural, suitable for tubes Ø 10 and 11 mm, 1,000 piece(s)/bag

**Product characteristics**

Cap	standard cap
Closure type	push cap
Suitable for	tubes Ø 10 and 11 mm

**Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	natural

**Purity & certification**

Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Push cap, natural, suitable for tubes Ø 10 and 11 mm**

**Page 2**

**Packaging**

Minimum order qty.	10000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	10000
Piece(s) / pallet	600000
Depth of case	398 mm
Width of case	358 mm
Height of case	175 mm
Case volume	0.0249 cbm
Weight of product	0.0004 kg
Weight of case	4.9 kg
EAN of inner box	4038917084093
EAN case	4038917002950

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**