



### Product description

Order number	65.793.516
Product description	Push cap, orange, suitable for tubes Ø 16-17 mm, 1,000 piece(s)/bag

### Product characteristics

Closure type	push cap
Suitable for	tubes Ø 16-17 mm

### Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	orange

### Purity & certification

Product category	no medical device   IVD
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	200000
Case volume	0.0414 cbm
Weight of product	0.001 kg
Weight of case	7.62 kg
EAN of inner box	4038917120869
EAN case	4038917033718

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

