



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 65.809.303   |
| Product description | Push cap, violet, suitable for tubes Ø 12 mm, 1,000 piece(s)/bag |

### **Product characteristics**

|              |               |
|--------------|---------------|
| Closure type | push cap      |
| Suitable for | tubes Ø 12 mm |

### **Material & colours**

|                   |                                  |
|-------------------|----------------------------------|
| Product material  | Low Density Polyethylene (LD-PE) |
| Colour of product | violet                           |

### **Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | yes                     |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 1000          |
| Piece (s) / outer case        | 10000         |
| Piece(s) / pallet             | 400000        |
| Case volume                   | 0.0382 cbm    |
| Weight of product             | 0.0006 kg     |
| Weight of case                | 6.82 kg       |
| EAN of inner box              | 4038917369053 |
| EAN case                      | 4038917369046 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

