



Product Specification
Push cap, natural, suitable for tubes Ø 11.5 and 12 mm

Page 1



Product description

Order number	65.809
Product description	Push cap, natural, suitable for tubes Ø 11.5 and 12 mm, 1,000 piece(s)/bag

Product characteristics

Cap	standard cap
Closure type	push cap
Suitable for	tubes Ø 11.5 and 12 mm

Size

Diameter	13.2 mm
Height of product	15.1 mm

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	natural

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Push cap, natural, suitable for tubes Ø 11.5 and 12 mm

Page 2

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	10000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	10000
Piece(s) / pallet	420000
Depth of case	398 mm
Width of case	398 mm
Height of case	250 mm
Case volume	0.0396 cbm
Weight of product	0.0006 kg
Weight of case	6.66 kg
EAN of inner box	4038917084239
EAN case	4038917003094

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-15

This document was prepared by EDP support and is valid without signature.



SARSTEDT