



Product description

| | |
|---------------------|--|
| Order number | 55.461 |
| Product description | Tube, working volume: 14 ml, (LxØ): 119 x 15.5 mm, material: PS, round base, transparent, push cap, 500 piece(s)/bag |

Product characteristics

| | |
|--------------------------|------------|
| Graduation | no |
| Closure type | push cap |
| Base shape | round base |
| Volume of work | 14 ml |
| Centrifugation max (RCF) | 3500 x g |

Size

| | |
|-------------------|---------|
| Diameter | 15.5 mm |
| Length of product | 119 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tube, 14 ml, (LxØ): 119 x 15.5 mm, PS

Page 2

Material & colours

| | |
|-------------------|------------------|
| Product material | Polystyrene (PS) |
| Colour of product | transparent |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 30000 |
| Depth of case | 292 mm |
| Width of case | 464 mm |
| Height of case | 413 mm |
| Case volume | 0.0559 cbm |
| Weight of product | 0.0058 kg |
| Weight of case | 6.29 kg |
| EAN of inner box | 4038917081696 |
| EAN case | 4038917000802 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-12

This document was prepared by EDP support and is valid without signature.



SARSTEDT