



Product description

| | |
|---------------------|---|
| Order number | 55.467.300 |
| Product description | Tube, working volume: 8 ml, (LxØ): 100 x 13 mm, material: PS, round base, transparent, push cap, 500 piece(s)/StackPack |

Product characteristics

| | |
|--------------------------|------------|
| Graduation | no |
| Closure type | push cap |
| Base shape | round base |
| Volume of work | 8 ml |
| Centrifugation max (RCF) | 3500 x g |

Size

| | |
|-------------------|--------|
| Diameter | 13 mm |
| Length of product | 100 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tube, 8 ml, (LxØ): 100 x 13 mm, PS

Page 2

Material & colours

| | |
|-------------------|------------------|
| Product material | Polystyrene (PS) |
| Colour of product | transparent |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | StackPack |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 2000 |
| Piece(s) / pallet | 96000 |
| Case volume | 0.0366 cbm |
| Weight of product | 0.0039 kg |
| Weight of case | 8.86 kg |
| EAN of inner box | 4038917367042 |
| EAN case | 4038917367035 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-12

This document was prepared by EDP support and is valid without signature.



SARSTEDT