



Product description

Order number	55.475.200
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PS, round base, transparent, push cap, suitable for Tecan Genesis, FE500, 490 piece(s)/StackPack

Product characteristics

Graduation	no
Closure type	push cap
Suitable for	Tecan Genesis, FE500
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	4000 x g

Size

Diameter	13 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tube, 5 ml, (LxØ): 75 x 13 mm, PS

Page 2

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1960
Type of smallest subpackaging	StackPack
Piece(s) / inner box	490
Piece (s) / outer case	1960
Piece(s) / pallet	109760
Depth of box	405 mm
Width of box	77 mm
Height of box	193 mm
Depth of case	428 mm
Width of case	383 mm
Height of case	175 mm
Case volume	0.0287 cbm
Weight of product	0.003 kg
Weight of case	6.83 kg
EAN of inner box	4038917158466
EAN case	4038917158459

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT