



### **Product description**

|                     |   |
|---------------------|---|
| Order number        | 55.476.013  |
| Product description | Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PS, round base, transparent, push cap, natural, sterile, 25 piece(s)/bag |

### **Product characteristics**

|                          |                 |
|--------------------------|-----------------|
| Graduation               | no              |
| Cap                      | Ventilation cap |
| Closure type             | push cap        |
| Base shape               | round base      |
| Volume of work           | 5 ml            |
| Centrifugation max (RCF) | 4000 x g        |

### **Size**

|                   |       |
|-------------------|-------|
| Diameter          | 12 mm |
| Length of product | 75 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tube, 5 ml, (LxØ): 75 x 12 mm, PS**

**Material & colours**

|                   |                                  |
|-------------------|----------------------------------|
| Product material  | Polystyrene (PS)                 |
| Colour of product | transparent                      |
| Closure material  | Low Density Polyethylene (LD-PE) |
| Colour of cap     | natural                          |

**Purity & certification**

|                           |                          |
|---------------------------|--------------------------|
| Satisfies the requirement | sealed at 0.6 bar/10 min |
| Product category          | no medical device   IVD  |
| Purity standard           | sterile                  |
| Sterilisation             | Electron irradiation     |
| Batched                   | yes                      |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 25            |
| Piece (s) / outer case        | 1000          |
| Piece(s) / pallet             | 48000         |
| Depth of case                 | 598 mm        |
| Width of case                 | 398 mm        |
| Height of case                | 145 mm        |
| Case volume                   | 0.0345 cbm    |
| Weight of product             | 0.0034 kg     |
| Weight of case                | 4.08 kg       |
| EAN of inner box              | 4038917107327 |
| EAN case                      | 4038917011785 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

