



Product description

Order number	55.484
Product description	Tube, working volume: 3.5 ml, (LxØ): 55 x 12 mm, material: PS, round base, transparent, push cap, 500 piece(s)/StackPack

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	3,5 ml
Centrifugation max (RCF)	5000 x g

Size

Diameter	12 mm
Length of product	55 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	4000
Type of smallest subpackaging	StackPack
Piece(s) / inner box	500
Piece (s) / outer case	4000
Piece(s) / pallet	168000
Depth of box	381 mm
Width of box	57 mm
Height of box	176 mm
Depth of case	388 mm
Width of case	378 mm
Height of case	250 mm
Case volume	0.0367 cbm
Weight of product	0.0019 kg
Weight of case	8.66 kg
EAN of inner box	4038917081917
EAN case	4038917001021

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

