



### **Product description**

Order number	55.522.300
Product description	Tube, working volume: 6.5 ml, (LxØ): 85 x 13 mm, material: PP, round base, transparent, push cap, 1,000 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	6,5 ml
Centrifugation max (RCF)	5000 x g
Autoclavable (max. °C)	121 °C

### **Size**

Diameter	13 mm
Length of product	85 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	72000
Case volume	0.0492 cbm
Weight of product	0.0029 kg
Weight of case	6.18 kg
EAN of inner box	4038917367226
EAN case	4038917367219

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

