



Product description

Order number	55.525.302
Product description	Carrier tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PP, round base, brown, push cap, 250 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	5000 x g
Autoclavable (max. °C)	121 °C

Size

Diameter	13 mm
Height of product	75 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Carrier tube, 5 ml, (LxØ): 75 x 13 mm, PP

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	brown

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	250
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	2000
Piece(s) / pallet	60000
Case volume	0.046 cbm
Weight of product	0.0026 kg
Weight of case	5.72 kg
EAN of inner box	4038917367257
EAN case	4038917367240

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT