




Product Specification
Tube, 5 ml, (LxØ): 75 x 12 mm, PP

Page 1



Product description

Order number	55.526.002
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PP, round base, brown, push cap, 500 piece(s)/StackPack

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	5000 x g
Autoclavable (max. °C)	121 °C

Size

Diameter	12 mm
Height of product	75 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tube, 5 ml, (LxØ): 75 x 12 mm, PP

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	brown

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	4000
Type of smallest subpackaging	StackPack
Piece(s) / inner box	500
Piece (s) / outer case	4000
Piece(s) / pallet	120000
Case volume	0.046 cbm
Weight of product	0.0021 kg
Weight of case	10.48 kg
EAN of inner box	4038917109284
EAN case	4038917031974

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT