




**Product Specification**  
**Tube, 5 ml, (LxØ): 75 x 12 mm, PP**

**Page 1**



**Product description**

Order number	55.526.005
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PP, round base, transparent, push cap, 1,000 piece(s)/bag

**Product characteristics**

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	5000 x g

**Size**

Diameter	12 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tube, 5 ml, (LxØ): 75 x 12 mm, PP**

**Page 2**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	4000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	4000
Piece(s) / pallet	80000
Depth of case	598 mm
Width of case	398 mm
Height of case	350 mm
Case volume	0.0833 cbm
Weight of product	0.0023 kg
Weight of case	10.15 kg
EAN of inner box	4038917109307
EAN case	4038917012478

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-21

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**