



### **Product description**

Order number	55.526.045
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PP, round base, transparent, push cap, with graduations, 2,000 piece(s)/bag

### **Product characteristics**

Graduation	yes
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	5000 x g
Autoclavable (max. °C)	121 °C

### **Size**

Diameter	12 mm
Height of product	75 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tube, 5 ml, (LxØ): 75 x 12 mm, PP**

**Page 2**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	4000
Type of smallest subpackaging	bag
Piece(s) / inner box	2000
Piece (s) / outer case	4000
Piece(s) / pallet	64000
Case volume	0.0813 cbm
Weight of product	0.002 kg
Weight of case	10.02 kg
EAN of inner box	4038917109635
EAN case	4038917068512

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**