



Product description

Order number	55.535
Product description	Tube, working volume: 3.5 ml, (LxØ): 55 x 12 mm, material: PP, round base, transparent, push cap, 1,000 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	3,5 ml
Centrifugation max (RCF)	5000 x g

Size

Diameter	12 mm
Length of product	55 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tube, 3.5 ml, (LxØ): 55 x 12 mm, PP

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	120000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0016 kg
Weight of case	8.96 kg
EAN of inner box	4038917082280
EAN case	4038917001380

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-08

This document was prepared by EDP support and is valid without signature.



SARSTEDT