



Product description

| | |
|---------------------|--|
| Order number | 57.462.015 |
| Product description | Tube, working volume: 12 ml, (LxØ): 110 x 17 mm, material: PS, conical base, transparent, push cap, with graduations, 250 piece(s)/StackPack |

Product characteristics

| | |
|--------------------------|--------------|
| Graduation | yes |
| Closure type | push cap |
| Base shape | conical base |
| Volume of work | 12 ml |
| Centrifugation max (RCF) | 2500 x g |

Size

| | |
|-------------------|--------|
| Diameter | 17 mm |
| Length of product | 110 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|------------------|
| Product material | Polystyrene (PS) |
| Colour of product | transparent |

Purity & certification

| | |
|---------------------------|--------------------------------|
| Satisfies the requirement | sealed at 0.6 bar/10 min |
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | StackPack |
| Piece(s) / inner box | 250 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 42000 |
| Depth of box | 366 mm |
| Width of box | 114 mm |
| Height of box | 201 mm |
| Depth of case | 388 mm |
| Width of case | 378 mm |
| Height of case | 250 mm |
| Case volume | 0.0367 cbm |
| Weight of product | 0.0063 kg |
| Weight of case | 6.96 kg |
| EAN of inner box | 4038917381079 |
| EAN case | 4038917381062 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

