



### Product description

Order number	58.490
Product description	Tube, working volume: 23 ml, (LxØ): 75 x 23.5 mm, material: PS, flat base, transparent, push cap, without cap, 500 piece(s)/bag

### Product characteristics

Graduation	no
Closure type	push cap
Base shape	flat base
Volume of work	23 ml
Centrifugation max (RCF)	1500 x g

### Size

Diameter	23.5 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tube, 23 ml, (LxØ): 75 x 23.5 mm, PS**

**Page 2**

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

**Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	24000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0065 kg
Weight of case	7.15 kg
EAN of inner box	4038917081979
EAN case	4038917001083

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-16

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**