



Product description

Order number	58.505.026
Product description	Tube, working volume: 12 ml, (LxØ): 95 x 16.5 mm, material: PP, flat base, transparent, push cap, natural, cap enclosed, 100 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	flat base
Volume of work	12 ml

Size

Diameter	16.5 mm
Length of product	95 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tube, 12 ml, (LxØ): 95 x 16.5 mm, PP

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Colour of cap	natural

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	28000
Case volume	0.0649 cbm
Weight of product	0.0066 kg
Weight of case	7.12 kg
EAN of inner box	4038917108522
EAN case	4038917032407

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT