



### **Product description**

Order number	58.536
Product description	Tube, working volume: 7 ml, (LxØ): 50 x 16 mm, material: PP, flat base, transparent, push cap, 1,000 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	push cap
Base shape	flat base
Volume of work	7 ml
Centrifugation max (RCF)	4000 x g

### **Size**

Diameter	16 mm
Length of product	50 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tube, 7 ml, (LxØ): 50 x 16 mm, PP**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	72000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0023 kg
Weight of case	5.29 kg
EAN of inner box	4038917082297
EAN case	4038917001397

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

