



Product Specification

Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP, without cap, 1,000 piece(s)/bag

Page 1



Product description

Order number	60.504.031
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 75 x 13 mm, round base, transparent, material: PP, without cap, 1,000 piece(s)/bag, 2,000 piece(s)/carton

Product characteristics

Graduation	no
Closure type	screw cap
Base shape	round base
Volume of work	5 ml

Size

Diameter	13 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP, without cap, 1,000 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	72000
Depth of case	458 mm
Width of case	408 mm
Height of case	260 mm
Case volume	0.0486 cbm
Weight of product	0.0026 kg
Weight of case	5.75 kg
EAN of inner box	4038917412735
EAN case	4038917412728

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

