



### Product description

Order number	60.541.009
Product description	Screw cap tube, working volume: 13 ml, (LxØ): 101 x 16.5 mm, material: PP, round base with skirted base, transparent, screw cap, without cap, 500 piece(s)/bag

### Product characteristics

Graduation	no
Closure type	screw cap
Base shape	round base with skirted base
Volume of work	13 ml
Centrifugation max (RCF)	4000 x g
Autoclavable (max. °C)	121 °C

### Size

Diameter	16.5 mm
Length of product	101 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Tube material	Polypropylene (PP)

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	32000
Case volume	0.0558 cbm
Weight of product	0.0067 kg
Weight of case	7.2 kg
EAN of inner box	4038917110563
EAN case	4038917032438

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

