



### Product description

Order number	60.541.685
Product description	Screw cap tube, working volume: 13 ml, (LxØ): 101 x 16.5 mm, material: PP, round base with skirted base, transparent, screw cap, natural, cap enclosed, 500 piece(s)/bag

### Product characteristics

Graduation	no
Closure type	screw cap
Base shape	round base with skirted base
Volume of work	13 ml
Centrifugation max (RCF)	4000 x g

### Size

Diameter	16.5 mm
Length of product	101 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

Purity & certification

Satisfies the requirement	IATA
Product category	no medical device   IVD
Batched	no

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	24000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0066 kg
Weight of case	7.29 kg
EAN of inner box	4038917111065
EAN case	4038917013147

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

