



Product description

| | |
|---------------------|---|
| Order number | 60.542.307 |
| Product description | Screw cap tube, working volume: 8 ml, (LxØ): 57 x 16.5 mm, material: PP, flat base, transparent, screw cap, natural, cap enclosed, 500 piece(s)/bag |

Product characteristics

| | |
|--------------------------|-----------|
| Graduation | no |
| Closure type | screw cap |
| Base shape | flat base |
| Volume of work | 8 ml |
| Centrifugation max (RCF) | 4000 x g |
| Autoclavable (max. °C) | 121 °C |

Size

| | |
|-------------------|---------|
| Diameter | 16.5 mm |
| Length of product | 57 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|--|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | Polypropylene (PP) |
| Colour of cap | natural |
| O-ring material | Ethylene propylene diene rubber (EPDM) |

Purity & certification

| | |
|---------------------------|-------------------------|
| Satisfies the requirement | IATA |
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 40000 |
| Case volume | 0.0382 cbm |
| Weight of product | 0.0043 kg |
| Weight of case | 4.81 kg |
| EAN of inner box | 4038917376044 |
| EAN case | 4038917376037 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

