



### Product description

Order number	60.546.002
Product description	Screw cap tube, working volume: 6 ml, (LxØ): 92 x 11.5 mm, material: PP, round base, transparent, screw cap, natural, cap assembled, sterile, 500 piece(s)/bag

### Product characteristics

Graduation	no
Closure type	screw cap
Base shape	round base
Volume of work	6 ml
Centrifugation max (RCF)	4000 x g

### Size

Diameter	11.5 mm
Length of product	92 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

**Purity & certification**

Satisfies the requirement	IATA
Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	56000
Depth of case	598 mm
Width of case	398 mm
Height of case	125 mm
Case volume	0.0298 cbm
Weight of product	0.0031 kg
Weight of case	3.76 kg
EAN of inner box	4038917111768
EAN case	4038917013376

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

