



### Product description

|                     |  |
|---------------------|--|
| Order number        | 60.549.001   |
| Product description | Screw cap tube, working volume: 3.5 ml, (LxØ): 66 x 11.5 mm, material: PP, skirted conical base, transparent, screw cap, natural, cap assembled, sterile, 100 piece(s)/bag |

### Product characteristics

|                          |                      |
|--------------------------|----------------------|
| Graduation               | no                   |
| Cap                      | standard screw cap   |
| Closure type             | screw cap            |
| Base shape               | skirted conical base |
| Volume of work           | 3,5 ml               |
| Centrifugation max (RCF) | 4000 x g             |

### Size

|                   |         |
|-------------------|---------|
| Diameter          | 11.5 mm |
| Length of product | 66 mm   |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

|                   |                                   |
|-------------------|-----------------------------------|
| Product material  | Polypropylene (PP)                |
| Colour of product | transparent                       |
| Closure material  | High Density Polyethylene (HD-PE) |
| Colour of cap     | natural                           |

### **Purity & certification**

|                           |                                |
|---------------------------|--------------------------------|
| Satisfies the requirement | IATA                           |
| Product category          | In vitro diagnostic device, CE |
| CE certified              | CE                             |
| Purity standard           | sterile                        |
| Sterilisation             | Electron irradiation           |
| Batched                   | yes                            |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 100           |
| Piece (s) / outer case        | 1000          |
| Piece(s) / pallet             | 56000         |
| Depth of case                 | 598 mm        |
| Width of case                 | 398 mm        |
| Height of case                | 125 mm        |
| Case volume                   | 0.0298 cbm    |
| Weight of product             | 0.0022 kg     |
| Weight of case                | 2.81 kg       |
| EAN of inner box              | 4038917112017 |
| EAN case                      | 4038917054928 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

