



### **Product description**

Order number	60.550.020
Product description	Screw cap tube, working volume: 7 ml, (LxØ): 82 x 13 mm, material: PP, flat base, transparent, screw cap, without cap, 1,000 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	screw cap
Base shape	flat base
Volume of work	7 ml
Centrifugation max (RCF)	4000 x g
Autoclavable (max. °C)	121 °C

### **Size**

Diameter	13 mm
Height of product	82 mm
Length of product	82 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	3000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	3000
Piece(s) / pallet	48000
Case volume	0.0813 cbm
Weight of product	0.0029 kg
Weight of case	10.25 kg
EAN of inner box	4038917112284
EAN case	4038917032537

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

