



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 60.550.034   |
| Product description | Screw cap tube, working volume: 7 ml, (LxØ): 82 x 13 mm, material: PP, round base, transparent, screw cap, natural, cap enclosed, 1,000 piece(s)/bag |

### **Product characteristics**

|                          |            |
|--------------------------|------------|
| Graduation               | no         |
| Closure type             | screw cap  |
| Base shape               | round base |
| Volume of work           | 7 ml       |
| Centrifugation max (RCF) | 4000 x g   |
| Autoclavable (max. °C)   | 121 °C     |

### **Size**

|                   |       |
|-------------------|-------|
| Diameter          | 13 mm |
| Length of product | 82 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

|                   |                                   |
|-------------------|-----------------------------------|
| Product material  | Polypropylene (PP)                |
| Colour of product | transparent                       |
| Closure material  | High Density Polyethylene (HD-PE) |
| Colour of cap     | natural                           |

**Purity & certification**

|                           |                         |
|---------------------------|-------------------------|
| Satisfies the requirement | IATA                    |
| Product category          | no medical device   IVD |
| Batched                   | yes                     |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 1000          |
| Piece (s) / outer case        | 1000          |
| Piece(s) / pallet             | 56000         |
| Case volume                   | 0.029 cbm     |
| Weight of product             | 0.0037 kg     |
| Weight of case                | 4.17 kg       |
| EAN of inner box              | 4038917112390 |
| EAN case                      | 4038917013659 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

