



### **Product description**

Order number	60.550.100
Product description	Screw cap tube, working volume: 7 ml, (LxØ): 82 x 13 mm, material: PP, round base, transparent, screw cap, natural, cap enclosed, 1,000 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	screw cap
Base shape	round base
Volume of work	7 ml
Centrifugation max (RCF)	4000 x g

### **Size**

Diameter	13 mm
Diameter of opening	11,55 mm
Length of product	82 mm
Length including cap	82 mm
Length excluding cap	80.9 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

**Purity & certification**

Satisfies the requirement	IATA
Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	1000
Piece(s) / pallet	48000
Depth of case	388 mm
Width of case	378 mm
Height of case	215 mm
Case volume	0.0315 cbm
Weight of product	0.0034 kg
Weight of case	3.86 kg
EAN of inner box	4038917112574
EAN case	4038917013741

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

