



### Product description

Order number	60.557.010
Product description	Screw cap tube, working volume: 4.5 ml, (LxØ): 75 x 12 mm, material: PP, round base, transparent, screw cap, without cap, 1,000 piece(s)/bag

### Product characteristics

Graduation	no
Closure type	screw cap
Base shape	round base
Volume of work	4,5 ml
Centrifugation max (RCF)	4000 x g
Autoclavable (max. °C)	121 °C

### Size

Diameter	12 mm
Height of product	75 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	80000
Case volume	0.0382 cbm
Weight of product	0.0021 kg
Weight of case	5.15 kg
EAN of inner box	4038917113366
EAN case	4038917056441

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

