



## Product Specification

**Screw cap tube, 2.5 ml, 2.5 ml, (LxØ): 75 x 13 mm,  
conical false bottom, rounded tube bottom, PP,  
without cap, 500 piece(s)/StackPack**

**Page 1**



### Product description

|                     |                                                                                                                                                                                                                                       |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number        | 60.614.015                                                                                                                                                                                                                            |
| Product description | Screw cap tube, working volume: 2.5 ml, nominal volume: 2.5 ml, (LxØ): 75 x 13 mm, conical false bottom, rounded tube bottom, transparent, material: PP, with graduations, without cap, 500 piece(s)/StackPack, 2,000 piece(s)/carton |

### Product characteristics

|                |                                           |
|----------------|-------------------------------------------|
| Graduation     | yes                                       |
| Base shape     | conical false bottom, rounded tube bottom |
| Volume of work | 2,5 ml                                    |

### Size

|                   |        |
|-------------------|--------|
| Sample volume     | 2.5 ml |
| Diameter          | 13 mm  |
| Height of product | 75 mm  |
| Length of product | 75 mm  |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Screw cap tube, 2.5 ml, 2.5 ml, (LxØ): 75 x 13 mm,  
conical false bottom, rounded tube bottom, PP,  
without cap, 500 piece(s)/StackPack**

**Page 2**

### Material & colours

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | transparent        |
| Tube material     | Polypropylene (PP) |

### Purity & certification

|                  |                               |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified     | CE                            |
| Batched          | yes                           |

### Packaging

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 2000          |
| Type of smallest subpackaging | StackPack     |
| Piece(s) / inner box          | 500           |
| Piece (s) / outer case        | 2000          |
| Piece(s) / pallet             | 112000        |
| Depth of box                  | 413 mm        |
| Width of box                  | 76 mm         |
| Height of box                 | 189 mm        |
| Depth of case                 | 428 mm        |
| Width of case                 | 383 mm        |
| Height of case                | 175 mm        |
| Case volume                   | 0.0287 cbm    |
| Weight of product             | 0.0024 kg     |
| Weight of case                | 5.57 kg       |
| EAN of inner box              | 4038917280044 |
| EAN case                      | 4038917280037 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**