



Product description

| | |
|---------------------|---|
| Order number | 60.9921.830 |
| Product description | Screw cap tube, working volume: 10 ml, (LxØ): 97 x 16 mm, material: PS, skirted conical base, transparent, screw cap, yellow, cap assembled, 100 piece(s)/bag |

Product characteristics

| | |
|--------------------------|----------------------|
| Graduation | no |
| Closure type | screw cap |
| Base shape | skirted conical base |
| Volume of work | 10 ml |
| Centrifugation max (RCF) | 4000 x g |

Size

| | |
|-------------------|-------|
| Diameter | 16 mm |
| Length of product | 97 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | Polystyrene (PS) |
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | yellow |

Purity & certification

| | |
|---------------------------|--------------------------------|
| Satisfies the requirement | ADR |
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 30000 |
| Depth of case | 608 mm |
| Width of case | 324 mm |
| Height of case | 264 mm |
| Case volume | 0.052 cbm |
| Weight of product | 0.0079 kg |
| Weight of case | 8.7 kg |
| EAN of inner box | 4038917151139 |
| EAN case | 4038917027717 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

