



### Product description

Order number	60.9922.243
Product description	Screw cap tube, working volume: 25 ml, (LxØ): 90 x 25 mm, material: PP, skirted conical base, transparent, screw cap, yellow, cap assembled, sterile, 500 piece(s)/bag

### Product characteristics

Graduation	no
Closure type	screw cap
Base shape	skirted conical base
Volume of work	25 ml
Centrifugation max (RCF)	4500 x g

### Size

Diameter	25 mm
Length of product	90 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow

### **Purity & certification**

Satisfies the requirement	IATA
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

### **Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	500
Piece(s) / pallet	15000
Depth of case	608 mm
Width of case	324 mm
Height of case	264 mm
Case volume	0.052 cbm
Weight of product	0.0082 kg
Weight of case	4.86 kg
EAN of inner box	4038917151528
EAN case	4038917027762

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

