



### **Product description**

Order number	62.467.012
Product description	Tube, working volume: 8 ml, (LxØ): 100 x 13 mm, material: PS, round base, transparent, push cap, with print, label/print: white, with graduations, 500 piece(s)/bag

### **Product characteristics**

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	push cap
Base shape	round base
Volume of work	8 ml
Centrifugation max (RCF)	3500 x g

### **Size**

Diameter	13 mm
Height of product	100 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	56000
Case volume	0.0649 cbm
Weight of product	0.0041 kg
Weight of case	9.2 kg
EAN of inner box	4038917106856
EAN case	4038917032599

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

