



Product description

| | |
|---------------------|--|
| Order number | 62.476.178 |
| Product description | Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PS, round base, transparent, push cap, with print, label/print: white, 1,000 piece(s)/bag |

Product characteristics

| | |
|--------------------------|------------|
| Label/ Print | with print |
| Colour of print/label | white |
| Graduation | no |
| Closure type | push cap |
| Base shape | round base |
| Volume of work | 5 ml |
| Centrifugation max (RCF) | 3500 x g |

Size

| | |
|-------------------|-------|
| Diameter | 12 mm |
| Height of product | 75 mm |
| Length of product | 75 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|------------------|
| Product material | Polystyrene (PS) |
| Colour of product | transparent |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1000 |
| Piece (s) / outer case | 2000 |
| Piece(s) / pallet | 80000 |
| Case volume | 0.0382 cbm |
| Weight of product | 0.0026 kg |
| Weight of case | 6.2 kg |
| EAN of inner box | 4038917107662 |
| EAN case | 4038917011846 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

