



### Product description

Order number	62.492
Product description	Tube, working volume: 10 ml, (LxØ): 95 x 16.8 mm, material: PS, round base, transparent, push cap, with print, label/print: white, with graduations, 500 piece(s)/bag

### Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	push cap
Base shape	round base
Volume of work	10 ml
Centrifugation max (RCF)	4000 x g

### Size

Diameter	16.8 mm
Length of product	95 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

### **Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	36000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0058 kg
Weight of case	6.46 kg
EAN of inner box	4038917081993
EAN case	4038917001106

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

