



**Product Specification**  
**Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP, cap enclosed, 1,000 piece(s)/bag**

**Page 1**



**Product description**

Order number	62.504.038
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 75 x 13 mm, round base, brown, material: PP, with print, label/print: black, with graduations, cap enclosed, brown, 1,000 piece(s)/bag, 1,000 piece(s)/carton

**Product characteristics**

Label/ Print	with print
Colour of print/label	black
Graduation	yes
Closure type	screw cap
Light protection	yes
Base shape	round base
Volume of work	5 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP, cap enclosed, 1,000 piece(s)/bag**

**Size**

Diameter	13 mm
Length of product	75 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	brown
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown
Tube material	Polypropylene (PP)

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	1000
Piece(s) / pallet	56000
Depth of case	568 mm
Width of case	433 mm
Height of case	120 mm
Case volume	0.0295 cbm
Weight of product	0.0038 kg
Weight of case	4.39 kg
EAN of inner box	4038917412971
EAN case	4038917412964

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

