



### Product description

Order number	62.526.028
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PP, round base, transparent, push cap, natural, cap assembled, with print, label/print: white, with graduations, sterile, 1 piece(s)/blister

### Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Cap	Ventilation cap
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	5000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Size**

Diameter	12 mm
Length of product	75 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	natural

**Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	18000
Depth of case	578 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0415 cbm
Weight of product	0.0029 kg
Weight of case	2.82 kg
EAN of inner box	4038917109475
EAN case	4038917012492

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

