




Product Specification
Screw cap tube, 4.5 ml, (LxØ): 75 x 12 mm, PP,
with print

Page 1



Product description

Order number	62.557.012
Product description	Screw cap tube, working volume: 4.5 ml, (LxØ): 75 x 12 mm, material: PP, round base, transparent, screw cap, natural, cap assembled, with print, label/print: white, with graduations, sterile, 100 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	screw cap
Base shape	round base
Volume of work	4,5 ml
Centrifugation max (RCF)	5000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

Diameter	12 mm
Length of product	75 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

Purity & certification

Satisfies the requirement	IATA
Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	56000
Case volume	0.029 cbm
Weight of product	0.0046 kg
Weight of case	5.15 kg
EAN of inner box	4038917225762
EAN case	4038917225755

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

