



**Product Specification**  
**Screw cap tube, 50 ml, (LxØ): 115 x 28 mm, PP,**  
**with print**

**Page 1**



**Product description**

Order number	62.559.010
Product description	Screw cap tube, working volume: 50 ml, (LxØ): 115 x 28 mm, material: PP, skirted conical base, transparent, screw cap, natural, cap assembled, with print, label/print: blue/white, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 25 piece(s)/bag

**Product characteristics**

Label/ Print	with print
Colour of print/label	blue/white
Graduation	yes
Closure type	screw cap
Base shape	skirted conical base
Volume of work	50 ml
Centrifugation max (RCF)	8000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 50 ml, (LxØ): 115 x 28 mm, PP,  
with print**

**Page 2**

**Size**

Diameter	28 mm
Length of product	115 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

**Purity & certification**

Satisfies the requirement	ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

**Packaging**

Minimum order qty.	300
Type of smallest subpackaging	bag
Piece(s) / inner box	25
Piece (s) / outer case	300
Piece(s) / pallet	7200
Depth of case	588 mm
Width of case	398 mm
Height of case	270 mm
Case volume	0.0632 cbm
Weight of product	0.0146 kg
Weight of case	5.07 kg
EAN of inner box	4038917113571
EAN case	4038917066020

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

