



Product Specification
**Screw cap tube, 50 ml, (LxØ): 115 x 28 mm, PP,
with print**

Page 1



Product description

| | |
|---------------------|--|
| Order number | 62.559.102 |
| Product description | Screw cap tube, working volume: 50 ml, (LxØ): 115 x 28 mm, material: PP, skirted conical base, transparent, screw cap, red, cap enclosed, with print, label/print: white/blue, with graduations, 25 piece(s)/bag |

Product characteristics

| | |
|--------------------------|----------------------|
| Label/ Print | with print |
| Colour of print/label | white/blue |
| Graduation | yes |
| Closure type | screw cap |
| Base shape | skirted conical base |
| Volume of work | 50 ml |
| Centrifugation max (RCF) | 8000 x g |
| Autoclavable (max. °C) | 121 °C |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|-------------------|--------|
| Diameter | 28 mm |
| Length of product | 115 mm |

Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | red |

Purity & certification

| | |
|---------------------------|-------------------------|
| Satisfies the requirement | IATA |
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 300 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 25 |
| Piece (s) / outer case | 300 |
| Piece(s) / pallet | 8400 |
| Case volume | 0.0649 cbm |
| Weight of product | 0.0268 kg |
| Weight of case | 8.71 kg |
| EAN of inner box | 4038917376433 |
| EAN case | 4038917376426 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

