



Product Specification
Screw cap tube, 50 ml, (LxØ): 115 x 28 mm, PP,
with print

Page 1



Product description

Order number	62.559.205
Product description	Screw cap tube, working volume: 50 ml, (LxØ): 115 x 28 mm, material: PP, skirted conical base, transparent, screw cap, red, cap assembled, with print, label/print: white/blue, with graduations, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 25 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white/blue
Graduation	yes
Closure type	screw cap
Base shape	skirted conical base
Volume of work	50 ml
Centrifugation max (RCF)	8000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Screw cap tube, 50 ml, (LxØ): 115 x 28 mm, PP,
with print**

Page 2

Size

Diameter	28 mm
Length of product	115 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	red

Purity & certification

Satisfies the requirement	IATA
Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	300
Type of smallest subpackaging	bag
Piece(s) / inner box	25
Piece (s) / outer case	300
Piece(s) / pallet	8400
Case volume	0.0649 cbm
Weight of product	0.0155 kg
Weight of case	5.19 kg
EAN of inner box	4038917376402
EAN case	4038917376396

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

