



## Product description

Order number	62.9924.276
Product description	Tube, working volume: 10 ml, (LxØ): 100 x 16 mm, material: PP, conical base, transparent, push cap, without cap, with print, with graduations, 100 piece(s)/bag

## Product characteristics

Label/ Print	with print
Graduation	yes
Closure type	push cap
Base shape	conical base
Volume of work	10 ml
Centrifugation max (RCF)	4500 x g

## Size

Diameter	16 mm
Length of product	100 mm

## Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

## Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

## Packaging

Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	36000
Depth of case	260 mm
Width of case	415 mm
Height of case	465 mm
Case volume	0.0502 cbm
Weight of case	5 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tube, 10 ml, (LxØ): 100 x 16 mm, PP, with print**

**Page 2**

EAN of inner box	4038917932783
EAN case	4038917932790

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-20

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**