



**Product Specification**  
**Tube, 30 ml, (LxØ): 80 x 28 mm, PS, with paper**  
**label**

**Page 1**



**Product description**

Order number	63.9922.232
Product description	Tube, working volume: 30 ml, (LxØ): 80 x 28 mm, material: PS, flat base, transparent, screw cap, yellow, cap assembled, with paper label, label/print: yellow, sterile, 1 piece(s)/blister

**Product characteristics**

Label/ Print	with paper label
Colour of print/label	yellow
Graduation	no
Cap	standard screw cap
Closure type	screw cap
Base shape	flat base
Volume of work	30 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tube, 30 ml, (LxØ): 80 x 28 mm, PS, with paper label**

**Page 2**

**Size**

Diameter	28 mm
Length of product	80 mm

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Tube material	Polystyrene (PS)

**Purity & certification**

Satisfies the requirement	sealed at 18"hg/10min
Latex-free	yes
Purity standard	sterile
Pyrogen-free	no
Safety data sheet	no

**Packaging**

Minimum order qty.	300
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	300
Piece(s) / pallet	9000
Depth of case	548 mm
Width of case	373 mm
Height of case	408 mm
Case volume	0.0834 cbm
Weight of case	4.8 kg
EAN of inner box	4038917151412
EAN case	4038917037808

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

