



## Product Specification

**Multi-purpose container, 420 ml, (LxØ): 150 x 70 mm, graduated, PP, transparent, with safety label**

**Page 1**



### Product description

Order number	75.9922.818
Product description	Multi-purpose container, max. working volume: 420 ml, (LxØ): 150 x 70 mm, transparent, cap: yellow, graduated, material: PP, screw cap, with safety label, cap assembled, with safety label, sterile, 40 piece(s)/bag

### Product characteristics

Label/ Print	with safety label
Graduation	yes
Closure type	screw cap
Base shape	flat base
Volume of work	420 ml

### Size

Sample volume	420 ml
Diameter	70 mm
Height of product	150 mm
Length of product	150 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Multi-purpose container, 420 ml, (LxØ): 150 x 70 mm, graduated, PP, transparent, with safety label**

**Page 2**

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Tube material	Polypropylene (PP)

### Purity & certification

Product category	no medical device   IVD
Latex-free	yes
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

### Packaging

Minimum order qty.	80
Type of smallest subpackaging	bag
Piece(s) / inner box	40
Piece (s) / outer case	80
Piece(s) / pallet	1280
Depth of case	548 mm
Width of case	373 mm
Height of case	408 mm
Case volume	0.0834 cbm
Weight of product	0.0656 kg
Weight of case	6.03 kg
EAN of inner box	4038917240048
EAN case	4038917350952

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-28

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**