




Product Specification
Transfer pipette, 3.5 ml, (LxW): 155 x 15 mm,
LD-PE, transparent

Page 1



Product description

Order number	86.1171.010
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 155 x 15 mm, material: LD-PE, transparent, droplet size: 35 µl - 55 µl, injected graduation up to 3 ml, sterile, 10 piece(s)/blister

Product characteristics

Graduation	yes
------------	-----

Size

Sample volume	3.5 ml
Width of product	15 mm
Length of product	155 mm
Drop size	35 µl - 55 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Transfer pipette, 3.5 ml, (LxW): 155 x 15 mm,
LD-PE, transparent

Page 2

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Minimum order qty.	1200
Type of smallest subpackaging	blister
Piece(s) / inner box	300
Piece (s) / outer case	1200
Piece(s) / pallet	33600
Depth of box	261 mm
Width of box	161 mm
Height of box	231 mm
Depth of case	533 mm
Width of case	333 mm
Height of case	250 mm
Case volume	0.0444 cbm
Weight of product	0.0015 kg
Weight of case	3.06 kg
EAN of inner box	4038917353144
EAN case	4038917353113

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT