



Product Specification
Transfer pipette, 3.5 ml, (LxW): 156 x 12.5 mm,
LD-PE, transparent

Page 1



Product description

Order number	86.1172.001
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 156 x 12.5 mm, material: LD-PE, transparent, droplet size: 30 µl - 45 µl, injected graduation up to 1 ml, sterile, 1 piece(s)/blister

Product characteristics

Graduation	yes
------------	-----

Size

Sample volume	3.5 ml
Width of product	12.5 mm
Length of product	156 mm
Drop size	30 µl - 45 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Minimum order qty.	840
Type of smallest subpackaging	blister
Piece(s) / inner box	42
Piece (s) / outer case	840
Piece(s) / pallet	30240
Depth of case	578 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0415 cbm
Weight of product	0.0013 kg
Weight of case	2.42 kg
EAN of inner box	4038917126281
EAN case	4038917020565

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

