



Product Specification
Transfer pipette, 3.5 ml, (LxW): 156 x 12.5 mm,
LD-PE, transparent

Page 1



Product description

Order number	86.1172.020
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 156 x 12.5 mm, material: LD-PE, transparent, droplet size: 30 µl - 45 µl, sterile, 20 piece(s)/bag

Product characteristics

Graduation	yes
------------	-----

Size

Sample volume	3.5 ml
Width of product	12.5 mm
Length of product	156 mm
Drop size	30 µl - 45 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Transfer pipette, 3.5 ml, (LxW): 156 x 12.5 mm,
LD-PE, transparent

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	4000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	4000
Piece(s) / pallet	64000
Case volume	0.1309 cbm
Weight of product	0.0017 kg
Weight of case	8.64 kg
EAN of inner box	4038917318310
EAN case	4038917318303

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

