




Product Specification
Transfer pipette, 3.5 ml, (LxW): 155 x 12.5 mm,
LD-PE, transparent

Page 1



Product description

| | |
|---------------------|---|
| Order number | 86.1173.300 |
| Product description | Transfer pipette, nominal volume: 3.5 ml, (LxW): 155 x 12.5 mm, material: LD-PE, transparent, droplet size: 30 µl - 45 µl |

Product characteristics

| | |
|------------|----|
| Graduation | no |
|------------|----|

Size

| | |
|-------------------|---------------|
| Sample volume | 3.5 ml |
| Width of product | 12.5 mm |
| Length of product | 155 mm |
| Drop size | 30 µl - 45 µl |

Material & colours

| | |
|-------------------|----------------------------------|
| Product material | Low Density Polyethylene (LD-PE) |
| Colour of product | transparent |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|------------------------|---------------|
| Minimum order qty. | 500 |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 5000 |
| Piece(s) / pallet | 56000 |
| Case volume | 0.0994 cbm |
| Weight of product | 0.0012 kg |
| Weight of case | 6.89 kg |
| EAN of inner box | 4038917378475 |
| EAN case | 4038917378468 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

